


Automatic Seed Counter



MSASCCPLUSWISS Automatic Seed Counter: Specification:

- **Model:** MSASCCPLUSWISS
- **Count Range:** 1 to 12 mm
- **Count Accuracy:** $\pm 2\%$
- **Count Speed:** ?1000 seeds/3 minutes (for Sesame)
- **Counting Range:** 1 to 99,999 seeds
- **Noise Level:** <70 dB
- **Power Supply:** AC 220/110V, 50/60Hz
- **Operation Environment:** -10°C to 50°C, <85% RH
- **External Size (W x D x H):** 285 x 210 x 176 mm
- **Packing Size (W x D x H):** 430 x 280 x 230 mm
- **Gross Weight:** 9.5 kg

