


**Centrifuge Machine, Angle Rotor: 8 × 10 ml, Maximum Speed: 4000 rpm**



**MSBKCTL4CSWISS Centrifuge Machine:  
Specification:**

- **Model:** MSBKCTL4CSWISS
- **Angle Rotor:** 8 × 10 ml
- **Maximum Speed:** 4000 rpm
- **Maximum RCF:** 1880 × g
- **Maximum Capacity:** 8 × 20 ml
- **Speed Precision:** ±10 rpm
- **Time Range:** 1 second to 99 minutes 59 seconds
- **Power Supply:** AC 100~240V, 50/60Hz
- **Power Consumption:** 250W
- **Noise Level:** 75 dB
- **External Size (L × W × H):** 310 × 270 × 210 mm
- **Package Size (L × W × H):** 370 × 310 × 250 mm

- **Net Weight / Gross Weight:** 7.5 kg / 8.5 kg

This centrifuge is ideal for small-scale laboratory use with precise speed and time control.