

Analytical Balance, 220g, 4 Digit**?Electronic Analytical Balance (External Calibration) – Model:
CB220G4DICSWISS****Features:**

- **Humanized Design:** A novel, ergonomic structure with a fully transparent five-sided windshield for easy viewing and precise sample handling.
- **High Precision:** Equipped with a new-generation electromagnetic balance weighing sensor, ensuring accuracy in all measurements.
- **Advanced Applications:** Includes automatic fault detection, four-point linear calibration, and overload protection for robust operation.
- **Clear LCD Display:** Ultra-large LCD screen for clear, comfortable viewing during operation.

Functions:

- Enhanced Operational Features: Includes multiple functions such as auto-counting, unit conversion (e.g., metric carat, gold ounce), stability, full-scale taring, and zero-memory for streamlined processes.
- Connectivity: Built-in RS232C output interface for direct connection to external devices like computers and printers.
- Convenient Hook: Includes a built-in lower weighing hook for applications requiring suspended weighing.
- Power Supply:
- Voltage: AC110-220V, 50-60Hz
- Power Cord: Standard AC power supply; does not operate on batteries.

Technical Specifications:

- Model: CB220G4DICSWISS
- Capacity: 0 - 220g
- Readability: 0.1mg
- Scale Size: 80mm
- Stable Time: 6-8 seconds
- Repeatability: 0.3mg
- Linearity: 0.3mg
- Net Weight: 6.8 kg
- Dimensions: 351 x 216 x 341 mm
- Work Space Height: 240mm
- Calibration: External Calibration

This analytical balance is ideal for laboratories requiring precision and ease of operation, with external calibration capabilities to maintain accuracy.