


## Horizontal Electrophoresis Tank (Small Size)



### Horizontal Electrophoresis Tank (Small Size)

- **Model:** MSBKHET01SWISS
- **Gel Trays:**
  - 75mm x 60mm
- **Buffer Volume:** 300ml
- **Comb Thickness:** 1mm & 1.5mm
- **Comb Throughput:**
  - 6, 11 samples
- **Number of Samples:** 6~22
- **Standard Accessories:**
  - Body tank with electrode x1
  - Cap x1
  - Gel tray x1
  - Gel box x1

- 1.0mm thickness 6/11 wells comb x2
- 1.5mm thickness 6/11 wells comb x2