


Horizontal Electrophoresis Tank



Horizontal Electrophoresis Tank

- **Model:** MSBKHETO2SWISS
- **Gel Trays:**
 - 120mm x 120mm
 - 120mm x 60mm
 - 60mm x 120mm
 - 60mm x 60mm
- **Buffer Volume:** 700ml
- **Comb Thickness:** 1mm & 1.5mm
- **Comb Throughput:**
 - 2, 3, 6, 8, 11, 13, 18, 25 samples
- **Number of Samples:** 2~50
- **Standard Accessories:**
 - Body tank with electrode x1

- Cap x1
- Gel trays x5
- Gel box x1
- 1.5mm thickness 2/3 wells comb x1
- 1.0mm thickness 6/13 wells comb x1
- 1.0mm thickness 8/18 wells comb x1
- 1.0mm thickness 11/25 wells comb x4