

Large-Scale Electronic Balance, 60Kg



Large-Scale Electronic Platform Balance, 60Kg

Introduction:

The large-scale electronic balance combines a sleek, modern appearance with advanced design for faster and more efficient weighing. Its unique metal base and mechanical structure ensure superior stability and enhanced anti-interference capabilities. Designed with portability in mind, it features AC/DC exchangeable power for flexible use in field operations. Additionally, custom specifications are available upon request to meet specific user needs.

Features:

Stable Design: The innovative structure provides excellent stability, while the backlit LCD screen ensures comfortable viewing for operators.

Flexible Power Options: Supports both AC and DC power supplies, equipped with a valve-regulated lead-acid battery for portability.

Versatile Functions: Offers multiple measurement unit conversions, counting, percentage calculations, and multi-point calibration.

Connectivity: Built-in RS232C output interface enables seamless integration with computers, printers, and other devices.

Technical Parameters:

Model: MSCF60KG1GSWISS

Capacity: 0–60,000g

Readability: 10g

Scale Size: 500×400mm

Power Supply: 100–240V, 50/60Hz

External Dimensions (W×D×H): 532×302×682mm

Net Weight: 12kg

Package Size (W×D×H): 872×467×227mm

Gross Weight: 15kg

This balance is ideal for industrial, laboratory, and field operations where precision, durability, and connectivity are essential.