

Linear/Reciprocating Shaker



Feature Highlights

Speed: 40-200 rpm for all models.

Digital Panel: Speed and time control.

Working Pan Size: 25x25 cm / 32x32 cm, with a maximum load of 5 kg.

Operation Modes: Continuous and timed operation ranging from 1 minute to 23 hours 59 minutes.

Motor: DC brushless motor, ensuring safety, long life, and maintenance-free operation. Silent and stable during operation.

Technical Data

Model	:	MSLS40PROSWISS
Motor Type	:	Brushless motor
Plate Size	:	320 x 320 mm
Turn	:	Linear and reciprocating
Operation Mode	:	Linear and reciprocating 20 mm
Max Load	:	5 kg
Speed Range	:	40 - 200 rpm
Timing Range	:	1 min - 23 h 59 min
Display	:	LED display
Shell Protection Level	:	IP21
Power	:	45 W
Voltage	:	100-240 V, 50/60 Hz
Size	:	330 x 435 x 140 mm