


Refrigerated Centrifuge, Angle Rotor 12x15 ml, Speed 6000 rpm**MSBKCTL6CRSWISS Refrigerated Centrifuge:
Specification:**

- **Model:** MSBKCTL6CRSWISS
- **Rotor:** Angle Rotor 12x10/15 ml
- **Max. Speed:** 6000 rpm
- **Speed Accuracy:** ± 10 rpm
- **Max. RCF (Relative Centrifugal Force):** 5290 \times g
- **Max. Volume:** 4 \times 500 ml
- **Refrigeration System:** CFC-free compressor, R404a refrigerant
- **Timing Range:** 1 second to 99 hours 59 seconds / Continuous / Short-span
- **Temperature Range:** -20 $^{\circ}$ to +40 $^{\circ}$
- **Temperature Accuracy:** $\pm 1^{\circ}$
- **Noise Level:** ≤ 65 dB
- **Power Supply:** AC 220V, 50/60Hz

- **Rated Power:** 1200W
- **External Dimensions:** 768x672x444 mm

This centrifuge is equipped with a high-capacity rotor and a refrigeration system, ideal for applications requiring precise temperature control and speed accuracy.