

Semi-Micro Analytical Balance



Semi-Micro Analytical Balance

Electromagnetic Force Recovery Weighing System: Ensures precise and reliable measurements.

Highly Integrated Load Cell: Enhances accuracy and stability.

Recipe Database: Capable of storing up to 99 recipes, each containing up to 20 ingredients.

Large Graphic Display: Facilitates easy and efficient operation.

Multilingual Support: Available in English, German, French, Italian, Spanish, and Portuguese.

GLP/ISO Compliance: Records saved weight values with date, time, and serial number for compliance.

Large Windshield Glass Cover: Features three sliding doors for easy access and operation.

Keyboard Plastic Protection: Protects the keyboard from damage and contaminants.

Specifications

Model: MSBMA1265SWISS

Capacity: 62g/120g

Resolution: 0.01mg/0.1mg

Minimum Weighing: 60mg

Stable Time: ?6s

Display: LCD (white backlight with black font)

Pan Size: ?80mm

Operating Temperature: 5-35?, 20-80% humidity

Repeatability: 0.03mg/0.1mg

Linearity: $\pm 0.08\text{mg}/\pm 0.2\text{mg}$

Draft Shield Size: 152mm x 160mm x 225mm

Calibration Weight: Internal calibration

Optional Accessory: Printer

Package Dimensions: 345mm x 215mm x 345mm

Gross Weight/Net Weight: 8.25kg/5.8kg

Languages: Italian, English, Portuguese, German, French, Spanish

Packing Size: 510mm x 400mm x 620mm (without pallet)

Power Supply: 110-230V, 50/60Hz, Output: DC 12V 1A