

Ultrasonic Bath /Ultrasonic Cleaner, 6L



Ultrasonic Cleaner

Introduction

The Ultrasonic Cleaner is an innovative, non-polluting, and multi-purpose cleaning device, delivering exceptional results without damaging the objects being cleaned. It is particularly effective for items with complex shapes and accelerates processes such as dissolution and emulsification. With unparalleled advantages over traditional physical or chemical cleaning methods, the Ultrasonic Cleaner finds applications in a wide array of industries:

- Service industry

- Electronics industry
- Laboratories
- Machinery manufacturing
- Cemented carbide industry
- Chemical industry

Features

1. Standard Accessories

- Soundproof cover (1 pc)
- Basket (1 pc)
- PVC pipe (1 pc)

2. Construction

- Made with durable 304 stainless steel for both the shell and inner cavity, ensuring anti-corrosion properties.

3. Key Advantages of the Eco Series

- LED touch screen for simple and intuitive operation
- Digital timer supporting up to 199 minutes
- Adjustable heating capability ranging from RT to 80?
- Includes a tray and a lid for added convenience
- Stainless steel 304 tank with a powder-coated finish, meeting antiseptic and

hygienic standards for laboratory environments

- Equipped with safety handles for portability
- Maintains ultrasonic efficiency at varying temperatures
- Features frequency overload and thermal protection for enhanced safety

Technical Parameters

Specification	Details
Model	MSVD06ESWISS
Capacity (L)	6.5L
Ultrasonic Frequency	Single frequency, 40kHz
Time Range	Up to 199 minutes
Heating Temperature	RT to 80°
Ultrasonic Power (W)	180
Heating Power (W)	500
Power Supply	AC220V, 50Hz
Tank Size (mm)	300 × 150 × 150
External Size (mm)	330 × 180 × 310
Packing Size (mm)	430 × 250 × 350
Gross Weight (kg)	6
Drainage Outlet	Yes

This Ultrasonic Cleaner is an ideal solution for a wide range of cleaning applications, providing unmatched efficiency, reliability, and safety across diverse environments.