


## Vertical Electrophoresis Tank



### **MSBKVET01SWISS Vertical Electrophoresis Tank: Specifications:**

- **Model:** MSBKVET01SWISS
- **Glass Tray Dimensions (W\*L):**
  - Thick type: 101\*83 mm
  - Thin type: 101\*73 mm
- **Gel Trays Dimensions (W\*L):**
  - Hand glue: 83\*73 mm
  - Precast glue: 86\*68 mm
- **Comb Thickness:**
  - Standard: 1.0 mm
  - Optional: 1.5 mm
- **Comb Throughput:**

- 10, 15 samples
- **Amount of Gel:**
  - 1~2 pieces
- **Number of Samples:**
  - 10~30 pieces
- **Buffer Volume:**
  - 700 ml
- **Standard Accessories:**
  - Body tank \*1
  - Cap \*1
  - Electrode core \*1 set
  - Gel casting base \*2
  - Gel frame \*2
  - Single glue baffle \*1
  - Gel shovel \*1
  - Glass plates \*5 pairs
  - 1.0 mm thickness 10-well comb \*2
  - 1.0 mm thickness 15-well comb \*2

This vertical electrophoresis tank is designed for flexibility in gel thickness and comb throughput, making it ideal for running between 10 to 30 samples with consistent and reliable performance.