12-Channel Adjustable Volume Micro Pipettes

Model **MSPIPET8SWISS** | Brand **Master Science** | Unit **Pc**



12 channel Micropipettes, PIPE-T12 Description · Pipettes cover a volume range from 0.5µL to 300µL. · Ergonomic design provides excellent operating experience. · Large display window allows for easy volume identification. · Easy calibration and maintenance. · 12 channel pipettes are appropriate for 96 well plates. · Dispensing head rotates for effortless pipetting convenience. · Individual piston and tip cone assemblies allowing easy repair and maintenance.

KEY HIGHLIGHTS

- Compatible with most universal tip brands.

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

12 CHANNEL PIPETTES ARE APPROPRIATE FOR 96 WELL PLATES.

- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material-made tip cone secures high sealing performance.
- **Compatible with most universal tip brands.**

Technical Specifications

PARAMETER	SPECIFICATION
Volume RangeIncrement | Test Volume | (Accuracy error) | (Precision error)	
%µL | % | µL	
0.5-10µL0.1µL | 10µL | 1.50 | 0.15 | 1.50 | 0.15	
5µL2.50 | 0.125 | 2.50 | 0.125	
1µL4.00 | 0.04 | 4.00 | 0.04	
5-50µL0.5µL | 50µL | 1.00 | 0.5 | 0.50 | 0.25	
25µL1.50 | 0.375 | 1.00 | 0.25	
5µL3.00 | 0.15 | 2.00 | 0.1	
10-100µL1µL | 100µL | 0.80 | 0.8 | 0.30 | 0.3	
50µL1.00 | 0.5 | 0.80 | 0.4	

PARAMETER	SPECIFICATION
10µL3.00 | 0.3 | 2.00 | 0.2	
50-300µL5µL | 300µL | 0.70 | 2.1 | 0.25 | 0.75	
150µL1.00 | 1.5 | 0.50 | 0.75	
50µL1.50 | 0.75 | 0.80 | 0.4	