-25℃ Freezer, 100L

Model **MSBDF25V100SWISS** | Brand **Master Science** | Unit **Pc**

-25℃ Freezer, 100L

-25℃ Freezer, 100LSpecification Model: MSBDF25V100SWISS Capacity: 100L Type: Vertical Temperature Range: -10℃ to -25℃ Temperature Accuracy: 0.1℃ Control System: Microprocessor control with a digital display Alarm System: High and low temperature alarm Sensor failure alarm Door open alarm Refrigeration Type: Direct refrigeration Refrigerant: CFC-free Condenser: Bundy tube

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

GENERAL

General	Capacity100L
General	TypeVertical
General	Temperature Range-10℃ to -25℃
General	Temperature Accuracy0.1℃
General	Control SystemMicroprocessor control with a digital display
General	Alarm System (Refrigeration Type)Direct refrigeration
General	Alarm System (Refrigerant)CFC-free
General	Alarm System (Condenser)Bundy tube
General	High and low temperature alarm
General	Sensor failure alarm
General	Door open alarm

ADDITIONAL INFORMATION

Additional Information	25℃ Freezer, 100L