
Centrifuge Machine, Angle Rotor: 8 × 10 ml, Maximum Speed: 4000 rpm

Brand **Master Science** | Unit **Pc**



MSBKCTL4CSWISS Centrifuge Machine:Specification: Model: MSBKCTL4CSWISS Angle Rotor: 8 × 10 ml Maximum Speed: 4000 rpm Maximum RCF: 1880 × g Maximum Capacity: 8 × 20 ml Speed Precision: ±10 rpm Time Range: 1 second to 99 minutes 59 seconds Power Supply: AC 100~240V, 50/60Hz Power Consumption: 250W Noise Level: ≤55 dB External Size (L × W × H): 310 × 270 × 210 mm Package Size (L × W × H): 370 × 310 × 250 mm Net Weight / Gross Weight: 7.5 kg / 8.5 kg This centrifuge is ideal...

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Angle Rotor	8 × 10 ml
Maximum Speed	4000 rpm
Maximum RCF	1880 × g
Maximum Capacity	8 × 20 ml
Speed Precision	±10 rpm
Time Range	1 second to 99 minutes 59 seconds
Power Supply	AC 100~240V, 50/60Hz
Power Consumption	250W
Noise Level	≤55 dB
External Size (L × W × H)	310 × 270 × 210 mm

Standard Accessories Included

- 370 × 310 × 250 mm
- Net Weight / Gross Weight: 7.5 kg / 8.5 kg
- This centrifuge is ideal for small-scale laboratory use with precise speed and time control.