Classic Single-Channel Adjustable Micro Pipette

Model **MSCLBQ211SWISS** | Brand **Master Science** | Unit **Pc**



Micro Pipette

The Classic Single-Channel Adjustable Micro Pipette offers exceptional reliability and accuracy, making it a trusted choice for laboratories involved in precision liquid handling tasks.

KEY HIGHLIGHTS

- Adjustable Volume Range
- High Accuracy and Precision
- Digital/Scale Display
- Durability

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

Adjustable Volume Range
Facilitates liquid handling for a wide array of experimental needs by allowing precise adjustment of aspirated and dispensed liquid volumes.

High Accuracy and Precision
Ideal for tasks requiring stringent volume control to meet experimental requirements.

Digital/Scale Display
Simplifies volume reading and setting, enhancing operational efficiency.

ENHANCES USER COMFORT AND RELIABILITY DURING PROLONGED USE.

Durability
Resistant to impact, UV exposure, and chemical corrosion, ensuring long-term reliability.

Technical Specifications

PARAMETER	SPECIFICATION
Available in various models and ranges to cater to different liquid transfer applications.	
10 ~ 100 Micro Liter	
20 ~ 200 Micro Liter	
100 ~ 1000 Micro Liter	
1000 ~ 5000 Micro Liter	