8-Channel Adjustable Volume Micro Pipettes

Model **MSPIPET8SWISS** | Brand **Master Science** | Unit **Pc**



Dispensing head rotates for effortless pipetting convenience.

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

GENERAL

General	Volume RangeIncrement | Test Volume | (Accuracy error) | (Precision error)
General	%µL | % | µL
General	0.5-10µL0.1µL | 10µL | 1.50 | 0.15 | 1.50 | 0.15
General	5µL2.50 | 0.125 | 2.50 | 0.125
General	1µL4.00 | 0.04 | 4.00 | 0.04
General	5-50µL0.5µL | 50µL | 1.00 | 0.5 | 0.50 | 0.25
General	25µL1.50 | 0.375 | 1.00 | 0.25
General	5µL3.00 | 0.15 | 2.00 | 0.1
General	10-100µL1µL | 100µL | 0.80 | 0.8 | 0.30 | 0.3
General	50µL1.00 | 0.5 | 0.80 | 0.4
General	10µL3.00 | 0.3 | 2.00 | 0.2
General	50-300µL5µL | 300µL | 0.70 | 2.1 | 0.25 | 0.75
General	150µL1.00 | 1.5 | 0.50 | 0.75
General	50µL1.50 | 0.75 | 0.80 | 0.4

ADDITIONAL INFORMATION

Additional Information	Pipettes cover a volume range from 0.5µL to 1000µL.
Additional Information	Ergonomic design provides excellent operating experience.

ADDITIONAL INFORMATION

Additional Information	Large display window allows for easy volume identification.
Additional Information	Easy calibration and maintenance.
Additional Information	8 channel pipettes are appropriate for 96 well plates.
Additional Information	Individual piston and tip cone assemblies allowing easy repair and maintenance.
Additional Information	Compound material-made tip cone secures high sealing performance.
Additional Information	Compatible with most universal tip brands.