Electrophoresis Power Supply

Brand **Master Science** | Unit **Pc**



This power supply is designed to offer precise control for electrophoresis applications with a wide voltage, current, and power output range, making it suitable for various laboratory needs. The inclusion of parallel outputs for up to 4 groups provides enhanced flexibility for multiple runs.

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
2-300V (Resolution	1V)
5-2000mA (Resolution	1mA)
1-300W (Resolution	1W)
220V, 50Hz	
4 groups	
1 minute to 99 hours 59 minutes, adjustable	