
Large Scale Electronic Platform Balance, 150Kg

Model **MSBE150000LSWISS** | Brand **Master Science** | Unit **Pc**



Large-Scale Electronic Balance

MSBE150000LSWISS Large Scale Electronic Platform Balance Specification: • Model: MSBE150000LSWISS • Capacity: 0 to 150,000 g • Readability: 1 g • Scale Size: 500 × 400 mm • Power Supply: 100~240V, 50/60Hz • External Size (W × D × H): 530 × 300 × 680 mm • Net Weight: 11 kg • Package Size (W × D × H): 870 × 465 × 225 mm • Gross Weight: 14 kg

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Capacity	0 to 150,000 g
Readability	1 g
Scale Size	500 × 400 mm
Power Supply	100~240V, 50/60Hz
External Size (W × D × H)	530 × 300 × 680 mm
Net Weight	11 kg

Standard Accessories Included

- 870 × 465 × 225 mm
- Gross Weight: 14 kg