
Horizontal Electrophoresis Tank (Small Size)

Brand **Master Science** | Unit **Pc**



Horizontal Electrophoresis Tank (Small Size) • Model: MSBKHET01SWISS • Gel Trays: o 75mm x 60mm • Buffer Volume: 300ml • Comb Thickness: 1mm & 1.5mm • Comb Throughput: o 6, 11 samples • Number of Samples: 6~22 • Standard Accessories: o Body tank with electrode x1 o Cap x1 o Gel tray x1 o Gel box x1 o 1.0mm thickness 6/11 wells comb x2 o 1.5mm thickness 6/11 wells comb x2

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Buffer Volume	300ml
Comb Thickness	1mm & 1.5mm
Number of Samples	6~22

Additional Information

GEL TRAYS
75mm x 60mm
COMB THROUGHPUT
6, 11 samples

Standard Accessories Included

- Body tank with electrode x1
- Cap x1
- Gel tray x1
- Gel box x1
- 1.0mm thickness 6/11 wells comb x2
- 1.5mm thickness 6/11 wells comb x2