
Horizontal Electrophoresis Tank

Brand **Master Science** | Unit **Pc**



Horizontal Electrophoresis Tank • Model: MSBKHETO2SWISS • Gel Trays: o 120mm x 120mm o 120mm x 60mm o 60mm x 120mm o 60mm x 60mm • Buffer Volume: 700ml • Comb Thickness: 1mm & 1.5mm • Comb Throughput: o 2, 3, 6, 8, 11, 13, 18, 25 samples • Number of Samples: 2~50 • Standard Accessories: o Body tank with electrode x1 o Cap x1 o Gel trays x5 o Gel box x1 o 1.5mm thickness 2/3 wells comb x1 o 1.0mm thickness 6/13 wells comb x1 o 1.0mm thickness 8/18 wells comb x1 o 1.0mm...

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Buffer Volume	700ml
Comb Thickness	1mm & 1.5mm
Number of Samples	2~50

Additional Information

GEL TRAYS
120mm x 120mm
120mm x 60mm
60mm x 120mm
60mm x 60mm

COMB THROUGHPUT
2, 3, 6, 8, 11, 13, 18, 25 samples

Standard Accessories Included

- Body tank with electrode x1
- Cap x1
- Gel trays x5
- Gel box x1
- 1.5mm thickness 2/3 wells comb x1
- 1.0mm thickness 6/13 wells comb x1
- 1.0mm thickness 8/18 wells comb x1
- 1.0mm thickness 11/25 wells comb x4