# Vortex Mixer

Model **MSVX4SWISS** | Brand **Master Science** | Unit **Pc**








The Vortex Mixer is a compact and high-performance laboratory tool designed to provide quick and uniform mixing for various liquid samples. With five preset speed buttons, a maintenance-free brushless motor, and adjustable vortex speeds, this mixer ensures high repeatability and reliability in every experiment.

### KEY HIGHLIGHTS

- Works quietly (<52 dB) without disrupting lab operations.
- Five Preset Buttons
- Stored Programs
- Brushless Motor

**Request a Quotation**
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

## Key Benefits

### LOW-NOISE PERFORMANCE

**Works quietly (<52 dB) without disrupting lab operations.**

### KEY FEATURES

**Five Preset Buttons**
Quick and easy speed selection for varied mixing applications.

**Stored Programs**
Supports multi-purpose use, increasing efficiency and workflow.

**Brushless Motor**
Maintenance-free and long-lasting for consistent performance.

## Technical Specifications

| PARAMETER | SPECIFICATION |
|---|---|
| Adjustable Vortex SpeedEnables controlled mixing from 200 to 3000 rpm. | |
| Quiet OperationLow sound emission (<52 dB) for a noise-free lab environment. | |
| Speed Range200 - 3000 rpm | |
| Maximum Vortex Speed3500 rpm | |
| Timer0s - 99h 59min | |
| Amplitude3mm | |
| Sound Emission<52 dB | |
| Power SupplyAC100-230V, 50/60Hz | |
| Power Consumption45W | |
| Dimensions (W×D×H)220 × 300 × 170 mm | |
| Net Weight/Gross Weight7.9 kg / 9.8 kg | |
| Shipping Dimensions (L×W×H)400 × 380 × 290 mm | |

## Applications

**Chemical & Biological Research Labs**

**Medical & Pharmaceutical Testing**

**Sample Preparation & Analysis**

Optimize your lab workflow with the precision and efficiency of the Vortex Mixer!

**Order Now to experience seamless and controlled mixing!**