Water Distiller, Water Control Type

Brand **Master Science** | Unit **Pc**



This machine is ideal for producing distilled water in laboratories within the medical and health industry, chemical industry, and scientific research units.

KEY HIGHLIGHTS

- Automatic Functionality
- Operational Indicators
- Electrode-Free Principle

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

EQUIPPED WITH WATER-DEFICIENT AND POWER-OFF ALARMS.

Automatic Functionality
Automatically replenishes water and reheats.

Operational Indicators
Includes power heating instructions and water level display.

Electrode-Free Principle
Prevents power-off issues due to electrode fouling from long-term use.

Technical Specifications

PARAMETER	SPECIFICATION
Water Output	5L/h
Power	5kW
Electricity	AC220V 50Hz
Overall Dimension (D*H)	310mm × 740mm
Net Weight/Gross Weight	7.1kg / 7.8kg